



Wuhan Greentek Pty. Ltd. Engineering drawing



Greentek standard 32ch Layout



Pinout chart

Processing requirements:

- 1、 Cap material: Elastic mesh cap;
- 2、 Electrode material: Sintered Silver / Silver chloride;
- 3、 Electrode Holder: Silicone rubber;
- 4、 Lead wire material: multi-strand Silver plated copper foil wire
+bulletproof fibre + medical TPU outsider;
- 5、 Lead wire length: 1.5m;
- 6、 Connector: DB37 (male);
- 7、 Drop lead with touch proof male connector for HEOL、 HEOR、 VEOL、 VEOR、 A1、 A2

Product: GT Cap ST - 32ch

Drawing No.: GREENTEK/TD-05

Version: V1.1