


# Wuhan Greentek Pty. Ltd. Engineering drawing



## Greentek 32ch Layout - Brain Products

### Processing requirements:

- 1、 Cap material: Elastic mesh cap;
- 2、 Electrode material: Sintered Silver / Silver chloride;
- 3、 Electrode Holder: Ultra-thin integrated holder
- 4、 Lead wire material: multi-strand Silver plated copper foil wire  
+bulletproof fibre + medical TPU outsider;
- 5、 Lead wire length: 1.5m;
- 6、 Connector: BP50PIN (male);
- 7、 Extra IO as floating electrode;

**Product: GT Cap Pro - BP 32ch**

Drawing No.:GREENTEK/TD-07  
Version: V1.1

**Wuhan Greentek Pty. Ltd.  
Engineering drawing**

<b>Channel Number</b>	<b>Name</b>	<b>PIN Number</b>
1	FP1	1
2	FP2	27
3	F3	28
4	F4	3
5	C3	4
6	C4	30
7	P3	31
8	P4	6
9	O1	7
10	O2	33
11	F7	34
12	F8	9
13	T7	10
14	T8	36
15	P7	37
16	P8	12
17	FZ	13
18	CZ	39
19	PZ	40
20	IO	15
21	FC1	16
22	FC2	42
23	CP1	43
24	CP2	18
25	FC5	19
26	FC6	45
27	CP5	46
28	CP6	21
29	FT9	22
30	FT10	48
31	TP9	49
32	TP10	24
N	GND	25
N	REF	38

**Pinout Chart**

**Product: GT Cap Pro - BP 32ch**

Drawing No.:GREENTEK/TD-08  
Version: V1.0